



In Vivo Buffer pH5.0

Description

In Vivo Buffer pH5.0

Product Details:

ichorbio's *in vivo* buffer pH5.0 is formulated and tested for *in vivo* applications. They are low in endotoxin and have been tested in animal models for toxicity. They can be used with all ichorbio *in vivo* antibodies and biosimilars.

Size:

ichorbio's buffer is only available in one standard size:

100ml

Please contact us for custom pricing.

Concentration:

1x

Formulation:

0.2 ?M filtered. 10mM Histidine, 10mM Aspartic Acid .PH 5.0±0.2. Contains no stabilizers or preservatives.

Endotoxin:

? 1.0 EU/mg as determined by the LAL method

Storage:

Store at 2-8°C for up to 12 months.

Use:

For research use only (RUO): it is not for diagnostic or therapeutic procedures.